
# Pediatric Reduced Energy Defibrillation Electrodes Instructions for Use

70-02071-01 A

Copyright © 2019 Cardiac Science Corporation. All rights reserved.

For use with the following Cardiac Science Automated External Defibrillators (AEDs) only: Powerheart AED G3 9300 series.

**Caution: Restricted sale.**

Federal law restricts this device to sale by or on the order of a physician or other licensed practitioner.

AT THE HEART OF SAVING  
LIVES<sup>®</sup>

## Indications for Use

Cardiac Science Pediatric Reduced Energy Defibrillation Electrodes are intended for use with only the following Cardiac Science Automated External Defibrillators (AEDs): Powerheart AED G3 9300 series. These electrodes, when used with the AED, provide reduced defibrillation energy and are intended for use only on children or infants up to 8 years old, or up to 55 lbs (25 kg). If the child appears older or larger, use Cardiac Science adult defibrillation electrodes. Do not delay treatment to determine the child's exact age or weight.

The AED is intended for use by personnel who have been trained in its operation. The user should be qualified by training in basic life support or other physician-authorized emergency medical response. The AED is indicated for emergency treatment of victims exhibiting symptoms of sudden cardiac arrest who are unresponsive and not breathing. Post-resuscitation, if the victim is breathing, the AED should be left attached to allow for acquisition and detection of the ECG rhythm. If a shockable ventricular tachyarrhythmia recurs, the AED will charge automatically and advise the operator to deliver therapy.

## Contraindications for Use

Do not use with Cardiac Science 9100 series AEDs, FirstSave 9200 AEDs, or non-Cardiac Science AEDs.

## Warnings and Cautions

**WARNING! Improper function.**

For use only with the following Cardiac Science AEDs: Powerheart AED G3 9300 series. These pads will not function properly with any other AEDs or defibrillators.

**WARNING! Shock hazard.**

During shock delivery, an AED applies high voltage to the patient. Do not touch the patient, pads or any other equipment connected to the patient when defibrillating.

**WARNING! Improper connection.**

DO NOT pre-connect Pediatric electrodes to an AED while in standby mode.

**WARNING! Do not reuse pads.**

Used pads may not adhere properly to patient. Improper pad adhesion may result in skin burns. Improper pad adhesion may result in improper AED performance. Used pads may cause patient-to-patient contamination.

**WARNING! No pacing function.**

Pads are not intended for use in pacing.

**Caution: Electrode degradation.**

DO NOT open the package until ready to use the pads.

**Caution: Equipment damage.**

Remove the Cardiac Science pediatric electrode pads before using other manufacturers' defibrillation equipment.

**Caution: Improper equipment performance.**

Using pads that are damaged or expired, or with lead wires damaged, may result in improper AED performance. Examine the pads before use; the package should be sealed and the expiration date should not have passed. Do not use pads if expired.

**Note:** Store electrodes at room temperature.

These symbols can be found on the packaging or on the electrodes:

Cardiac Science pediatric defibrillation electrodes are disposable. Single patient use only. Do not reuse.

Dangerous voltage

Latex free

Date of manufacture

Attention required. Consult the accompanying documents or AED Operation and Service Manual for additional information

Use pads by this date.

Consult instructions for use

Lot number

Operating temperature

Waste Electronic Electrical Equipment (WEEE). Separate collection for waste electrical and electronic equipment.

Model number

Manufacturer

Authorized representative in the European Community

Mark of conformity

## Instructions for Use

**WARNING! AED not rescue ready.**

If the rescuer performs the steps below in a different sequence than outlined, the AED may give a “check electrodes” or “check pads” message. Unlike Cardiac Science adult defibrillation electrodes, the construction of our pediatric electrodes does not allow a Cardiac Science AED to check that there is an adequate connection and that the electrodes are ready to use prior to their being applied to the patient.

1. Remove clothing from the patient’s chest and make sure the patient’s chest is clean and dry.
2. Tear open the pediatric pads package and remove the pads.
3. Carefully peel one pad away from the liner.

**WARNING! Equipment damage.**


DO NOT pull on the lead wires to separate the pads from the blue liner.

**WARNING! Improper therapy delivery.**

Failure to remove blue liner could impact therapy delivery.

4. Standard pads placement (Figure 1).
  - a. Place one pad with the sticky side to the patient’s bare skin, on the upper right chest as shown in Figure 1. Press the pad firmly to the patient’s skin.
  - b. Carefully peel the second pad away from the liner.
  - c. Place the second pad on the patient’s lower left chest as shown in Figure 1.
5. Open the AED lid. If adult electrode pads are attached to the AED, disconnect them before connecting the pediatric pads.
6. Connect the pediatric pads to the AED. To attach, match the connector colors (red to red) from the pads to the AED and slide together until the connectors are firmly attached.
7. The AED begins patient analysis. Refer to the AED Operator’s and Service Manual for additional operating information.


**Figure 1: Standard pads placement (recommended)**

**Note:** Pads can be placed in either position on the patient.

8. Alternative pads placement (Figure 2)
  - a. Place one pad with the sticky side to the patient's bare skin, in the center of the chest, as shown in Figure 2. Press the pad firmly to the patient's skin.
  - b. Carefully peel the second pad away from the liner.
  - c. Place the second pad on the patient's back as shown in Figure 2.


**Figure 2: Alternative pads placement**

**Note:** Pediatric electrode pads are for short-term use only. One set of Cardiac Science Pediatric Electrode Pads can be used to monitor a patient for 2 hours and deliver up to 15 shocks.

To remove electrode pads, slowly peel back the edge of the pads away from the patient's skin. Dispose of the used pads.

### Attenuated Energy Values With Cardiac Science Pediatric Electrodes and STAR Biphasic Waveform

Typical waveform:  
Low energy (200 VE)  
50 ohm patient impedance

**Initial shock: Ultra Low Energy (150 VE) – All values are typical**

Patient's Impedance (ohms)	Phase 1		Phase 2		
	Voltage (volts)	Duration (ms)	Voltage (volts)	Duration (ms)	Energy (Joules)
25	370	6.1	258	3.2	31
50	550	7.3	366	3.2	36
75	640	8.6	417	3.2	37
100	705	9.8	442	3.2	36
125	770	11.1	453	3.2	35

**Initial shock: Low Energy (200 VE) – All values are typical**

Patient's Impedance (ohms)	Phase 1		Phase 2		
	Voltage (volts)	Duration (ms)	Voltage (volts)	Duration (ms)	Energy (Joules)
25	430	6.1	298	3.2	42
50	630	7.3	422	3.2	50
75	745	8.6	482	3.2	51
100	790	9.8	511	3.2	49
125	855	11.1	524	3.2	47

**Second and subsequent shocks: Ultra Low Energy (150 VE) – All values are typical**

Patient's Impedance (ohms)	Phase 1		Phase 2		
	Voltage (volts)	Duration (ms)	Voltage (volts)	Duration (ms)	Energy (Joules)
25	370	5.8	270	3.2	31
50	550	6.5	390	3.2	35
75	640	7.0	470	3.2	34
100	705	7.4	510	3.2	32
125	770	7.8	545	3.2	29

**Second and subsequent shocks: Low Energy (200 VE) – All values are typical**

Patient's Impedance (ohms)	Phase 1		Phase 2		
	Voltage (volts)	Duration (ms)	Voltage (volts)	Duration (ms)	Energy (Joules)
25	430	5.8	295	3.2	41
50	630	6.5	425	3.2	47
75	745	7.0	510	3.2	46
100	790	7.4	560	3.2	43
125	855	7.8	610	3.2	39

**Second and subsequent shocks: High Energy (300 VE) – All values are typical**

Patient's Impedance (ohms)	Phase 1		Phase 2		
	Voltage (volts)	Duration (ms)	Voltage (volts)	Duration (ms)	Energy (Joules)
25	500	5.8	380	3.2	56
50	700	6.5	520	3.2	63
75	820	7.0	620	3.2	62
100	920	7.4	680	3.2	58
125	960	7.8	720	3.2	53

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Copyright © 2019 Cardiac Science Corporation. Tous droits réservés.
Copyright © 2019 Cardiac Science Corporation. Se reservan todos los derechos.
Copyright © 2019 Cardiac Science Corporation. Todos os direitos reservados.
Copyright © 2019 Cardiac Science Corporation. Tutti i diritti riservati.
Copyright © 2019 Cardiac Science Corporation. Alle Rechte vorbehalten.
Copyright © 2019 Cardiac Science Corporation. Alle rechten voorbehouden.