

X Series®

Technical Specifications:



Because Medics Carry More These Days

- At less than 12 pounds/6 kilograms, the X Series is up to half the size and weight of other full-featured devices.
- Its large, bright screen allows for simultaneous viewing of four traces, all 12 ECG leads, or static and dynamic 12-lead views on a “split screen.”
- The latest, most advanced and most complete monitoring capabilities, including gold standard Masimo® full rainbow® SET and Welch Allyn’s Rapid NIBP technologies can be added.
- Trauma-specific parameters deliver insight on occult hemorrhage (SpHb™), fluid responsiveness (PVI™), and oxygenation status (SpOC™).

Because CPR Is a Focus for Cardiac Arrest

- CPR Dashboard™ provides real-time numeric depth and rate values and has a CPR timer that switches to an idle timer when compressions are stopped. The release indicator and the perfusion performance indicator (PPI) also help maximize CPR performance.
- See-Thru CPR® (unique to ZOLL®) helps minimize the duration of interruptions by filtering out artifact, allowing you to see organized underlying rhythms during CPR.

Simplify Charting with Advanced Communication

- The X Series is the first monitor/defibrillator with integrated WiFi. Bluetooth® and USB cellular modem capabilities are standard as well.
- With the push of a button, medics can now directly send X Series files to your system. With RescueNet CaseReview, access to case files is as fast as opening your web browser.
- The X Series acquires, stores, and transmits 12-leads to ZOLL’s easy-to-use and economical RescueNet 12-Lead and RescueNet MedGate 12-Lead management systems.

General

Weight: 11.7 lbs. (5.3 kg) with battery and paper

Dimensions: 8.9 in (22.6 cm) x 10.4 in (26.4 cm) x 7.9 in (20.1 cm)

Enclosure Protection: Solid Foreign Object: IEC 60529, IP5X

Water: IEC 60529, IPX5

Operating Temperature: 0 to 50°C

Humidity: 15 to 95% RH

(non-condensing)

Vibration: MIL-STD 810G, Method 514.6, 4.4.2 Procedure II

EN ISO 9919 (per IEC 60068-2.64)

EN 1789 for ambulance. RTCA/DO-160G (multiple helicopter frequencies)

Shock: IEC 60068-2-27, 100g, 6 ms half sine

Drop: EN 1789, 30-inch functional drop, IEC 60601-1, tested at 2 meters

Altitude: -170 m to 4572 m (-557 feet to 15,000 feet)

Monitor/Display

Size: 6.5 in (16.56 cm) diagonal

Input: 3-lead, 5-lead, or 12-lead patient cable, paddles, or multifunction electrodes

Type: Color LCD, 640 x 480 pixels, 800 MCD

Sweep Speed: 25 mm/sec or 50 mm/sec (user selectable)

Lead Selections: Paddles (Pads), I, II, III, AVR, AVL, AVF, V1-6

Frequency Response:

(user selectable)

0.67 – 20 Hz Limited mode

0.67 – 40 Hz Monitor mode

0.25 – 40 Hz Filtered diagnostic mode

0.05 – 150 Hz Diagnostic mode

Defibrillator

Waveform: Rectilinear Biphasic™

Energy Selections: 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 15, 20, 30, 50, 70, 85, 100, 120, 150, 200 joules

Energy Display: Display indicates selected and delivered energy

Charge Time: Less than 7 seconds with a new, fully charged battery

Charge Controls: Front panel and on apex paddle

Pads/Paddles: Hands-free pads and external adult and pediatric paddles. Adult paddles slide off to expose pediatric paddles.

Synchronized Mode: Synchronizes defibrillator pulse to patient's R wave.

"SYNC" message displayed on monitor and markers shown on both monitor and recorded ECG. Meets Clause 104 of AAMI DF80.

Patient Modes

User Selectable: Adult, Pediatric, Neonate.

Automatically sets configurable defaults for alarm limits, defibrillation energy, and NIBP settings.

Trends (on-screen)

Tabular numeric format

All parameters trended/viewable

Trend Intervals: 1, 5, 10, 15, 30, 60 minutes

Duration: 24 hours at 1-minute intervals

Snapshots: Minimum of 32 (12 seconds of pre- and post-button press)

Quick-access trend display via quick-access button

Memory Capacity

A combination of 24 hours of trends at 1-minute intervals, 1000 time-stamped events, and 32 snapshots, including monitor, defibrillator, pacer, and treatment snapshots

Pacer

Type: External transcutaneous pacing

Pulse Types: Rectilinear, constant current

Pulse Width: 40 ms ± 2 ms

Pacer Rate: 30 – 180 bpm ± 1.5%

Output Current: 0 – 140 mA ± 5% or 5 mA, whichever is greater

Modes: Demand and Fixed

CPR

CPR Feedback Technology

See-Thru CPR®: Artifact filtering

Real CPR Help®: Real-time depth and rate feedback

CPR Dashboard™: Numeric readout of depth and rate, release indicator, perfusion performance indicator (PPI)

ECG

Cable Detection: Automatic 3-, 5-, 12-lead ECG

Input: 3-lead cable, 5-lead cable,

12-lead cable, paddles, therapy cable

Leads: I, II, III, AVR, AVL, AVF, V1 – V6

Heart Rate Range: 30 – 300 bpm

Heart Rate Accuracy: ± 3 bpm or 3%

Pacer Detection and Display

ECG Sizes: 0.125, 0.25, 0.5, 1, 2, 4 cm/mV, and auto-ranging

Sweep Speed: 25 mm/sec, 50 mm/sec

Masimo® SET SpO₂

Saturation Range: 1 – 100%

Saturation Accuracy:

Oxygen saturation during no-motion conditions:

Adults, Pediatrics: 70 – 100%, ± 2 digits

0 – 69%, unspecified

Neonates: 70 – 100%, ± 3 digits

0 – 69%, unspecified

Oxygen saturation during motion conditions

Adults, Pediatrics: 70 – 100%, ± 3 digits

0 – 69%, unspecified

Neonates: 70 – 100%, ± 3 digits

0 – 69%, unspecified

Oxygen saturation during low-perfusion conditions:

Adults, Pediatrics: 70 – 100%, ± 2 digits

0 – 69%, unspecified

Neonates: 70 – 100%, ± 3 digits

0 – 69%, unspecified

Pulse Rate Range: 25 – 240 bpm

Pulse rate during no-motion conditions:

Adults, Pediatrics,

Neonates: 25 – 240 ± 3 digits

Pulse rate during motion conditions

Adults, Pediatrics,

Neonates: 25 – 240 ± 5 digits

SpO₂ Average time setting:

4, 8 (default), 16 seconds

Masimo rainbow® SET SpCO®

Range: 0 – 99%

Accuracy: 1 – 40% ± 3 digits

Masimo rainbow® SET SpMet®

Range: 0 – 99%

Accuracy: 1 – 15% ± 1 digits

Masimo rainbow® SET SpHb®

Measurement Range: 0 – 25 g/dL

Accuracy (Adults/Infants/Pediatrics): 8 – 17 g/dL +/- 1 g/dL

Resolution (SpHb g/dL): 0.1 g/dL

Masimo rainbow® SET SpOC™

Measurement Range: 0 – 35 ml of O₂/dL of blood

Resolution: 0.1 ml/dL

Masimo Perfusion Index (PI)

Measurement Range: 0.02 – 20%

Resolution: 0.1%

Masimo Pleth Variability Index (PVI®)

Measurement Range: 0 – 100%

Resolution: 1%

Non-Invasive Blood Pressure (NIBP)

Smartcuf® and Sure BP® NIBP technology

Measurement Intervals: Automatic 1-, 2-, 3-, 5-, 10-, 15-, 30-, 60-minute, and manual

quick-action NIBP Start/Stop button

TurboCuf: 5 min of repeated NIBP readings

Display: Systolic, diastolic, mean. Viewable

on-screen with large numerics.

Cuff Overpressure Protection

Typical Measurement Time:

30 – 45 sec (on deflation)

15 – 30 sec – Sure BP (on inflation)

Standard Cuff Sizes: Neonate #1 – #5, Infant, Small Child, Small Adult, Adult Long, Large Adult, Thigh

Default Cuff Inflation Pressure:

Adults: 160 mmHg

Pediatrics: 120 mmHg

Neonates: 90 mmHg

Pressure Measurement Range:

Systolic: 20 – 260 mmHg

Diastolic: 10 – 220 mmHg

Mean: 13 – 230 mmHg

Maximum Cuff Inflation Pressure:

Adults: 270 mmHg

Pediatrics: 170 mmHg

Neonates: 130 mmHg

Microstream® etCO₂

Range:

0 – 150 mmHg

Accuracy:

0 – 38: ± 2 mmHg

39 – 150: ± 5% + .08% mmHg >38

Respiration Range: 0 – 149 breaths per minute

Respiration Rate Accuracy:

0 – 70 bpm ± 1 bpm

71 – 120 bpm ± 2 bpm

121 – 149 bpm ± 3 bpm

Flow rate: 50 ml/min -7.5 + 15 ml/min

Typical response time: 2.9 seconds

Maximum response time: 3.9 seconds

Impedance Pneumography

Displayed Data: Numeric breath rate,

Impedance waveform

Breath Rate Range: Adult, Pediatric:

2 to 150 breaths/minute. Neonates: 3 to 150

breaths/minute.

Breath Rate Accuracy: 2% or ± 2, whichever is greater

Displayed Breath Rate: Average of last 10 breath-to-breath rates

Leads: Lead I (RA – LA), Lead II (RA – LL)

Sweep Speed: 3.13, 6.25, 12.5 mm/sec

Alarm Settings: High, low, and no breath rate alarm

Invasive Pressure

Three Channels

Pressure Range: -30 – 300 mmHg

Pulse Rate Measurement Range:

25 – 250 bpm

Formats: S/D, S/D (M), (M) user selectable

User-Selectable Labels: P1, P2, P3, ABP, AO, ART, CVP, BAP, FAP, LAP, PAP, RAP, UAP, UVP, ICP

Transducer Requirements:

5M μV/V/mmHg

Zero Adjustment: ± 200 mmHg

Numeric Pressure Accuracy: ± 2 mmHg or 2% reading, whichever is greater of reading, plus transducer error

Transducer Connector: Standard 6-pin MS connector

Temperature

Two YSI 400/700 series-compatible channels

Range: 0 – 50°C (32 – 122°F)

Units: °C or °F

Display: T1, T2, and Delta temp

Accuracy: ± 0.1°C (± 0.06°F) from 10 – 50°C (50 – 122°F) ± 0.2°C (± 0.11°F) from 0 – 10°C (32 – 50°F)

Printer

Type: High-resolution thermal array

Annotation: Time, date, ECG lead, ECG gain, heart rate, defibrillation, and pacing parameters and treatment summary events

Paper Width: 80 mm

Paper Speed: 25 mm/sec, 50 mm/sec

Delay: 6 seconds

Frequency Response: Automatically set to monitor's frequency response

Record Modes: Manual and automatic (user configurable)

Print Option Single waveform or a combination up to 3, on alarm, snapshots, treatment summary report, and trend summary

Battery

Type: Rechargeable lithium-ion, 11.1 Vdc, 6.6 Ah, 73 Wh

Capacity: With a new, fully charged battery operating at room temperature: At least 6 hours of continuous monitoring of ECG, SpO₂, CO₂, 3 invasive pressure channels, and 2 channels of temperature, with NIBP measurements every 15 minutes and 10 200-joule shocks. 300 shocks at 200 joules with fully charged new battery at room temperature. At least 3.5 hours pacing, with ECG, SpO₂, CO₂, three invasive pressures, temperature, NIBP every 15 minutes and pacing at 180 ppm, and 140 mA.

AC Power Adapter: 100 – 240 VAC 50, 60 Hz, 2A 100 – 115 VAC 400Hz, 2A

Specifications subject to change without notice.

© 2015 ZOLL Medical Corporation. "Advancing Resuscitation. Today.", CPR Dashboard, Rectilinear Biphasic, Real CPR Help, RescueNet, See-Thru CPR, X Series, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. Smartcuf, and Sure BP are registered trademarks of Welch Allyn. Masimo, Rainbow, SET, SpCO, SpHb, SpMet, and SpOC are trademarks or registered trademarks of Masimo Corporation. Microstream is a registered trademark of Oridion Medical 1987 Ltd. All other trademarks are the property of their respective owners.

Printed in U.S.A. MCN EP 1508 0108

**ADVANCING
RESUSCITATION.
TODAY.®**

ZOLL Medical Corporation Worldwide Headquarters

269 Mill Road
Chelmsford, MA 01824
978-421-9655
800-804-4356

For subsidiary addresses and fax numbers, as well as other global locations, please go to www.zoll.com/contacts.

ZOLL®