

Mobilize™

Compact bleeding control and airway management system



Technical Specifications

Be Ready for the Unexpected

When medical emergencies occur, the faster care can be provided, the better the outcome. Bystanders are frequently a victim's best hope for survival, but they don't always have the training or equipment they need to help. The Mobilize Compact bleeding control and airway management system includes medical supplies and a simple instructional app so even untrained bystanders can help provide care in the critical minutes before EMS arrives.

Modern Medical Supplies

Traditional first aid kits can help patch up scraped knees, but they don't generally include lifesaving supplies like tourniquets, pressure dressings, or face shields. The Compact bleeding control and airway management system is designed to address life-threatening emergencies such as severe bleeding, airway obstructions, and other traumatic injuries that can become fatal without intervention.

Real-time Instructions at Your Fingertips

The Mobilize Compact bleeding control and airway management system pairs with the Mobilize Rescue app, which users download and access from their smartphone. The app provides real-time guidance through simple interactive prompts, written instructions, and visual cues that show rescuers how to use the medical supplies in the kit to provide care during a medical emergency. The kit includes a multi-user app license so anyone with the Mobilize Rescue app will be prepared to help.

Have Peace of Mind at Work, at Home, or on the Road

Compact and lightweight, this bleeding control and airway management kit is handy to keep at work or take with you on the go. Military-grade nylon casing protects supplies in a range of environments and Molle straps allow you to easily attach it to a rucksack, belt, or other travel pack. With the Compact bleeding control and airway management system, you can have confidence that you — and anyone with the Mobilize Rescue app — are prepared to help during the average 7–14 minutes it takes for an ambulance to arrive.¹

ZOLL®

System Contents

Trauma Supplies

SOF-T Wide Tourniquet	1
QuikClot Bleeding Control Dressing	1
Hyfin Chest Seal	1
4" Emergency Trauma Dressing	1
CPR Face Shield With Bite Block	1
Emergency Space Mylar Blanket	1
Trauma Shears	1
Nitrile Gloves	8
Permanent Marker	1

Mobilize Rescue app with multi-user app license.

Treatment Capabilities

SEVERE BLEEDING

UNCONSCIOUSNESS

TROUBLE BREATHING

OPIATE OVERDOSES

HYPOTHERMIA

SEIZURES

CHOKING

CPR INSTRUCTION

AED APPLICATION

Mobilize™ Compact bleeding control and airway management system

Equipment

Modern medical supplies and a real-time instructional application to manage severe bleeding, airway complications, cardiac arrest, seizures, choking, and more.

Case Features


Military-grade nylon casing and a small, lightweight design make this kit both durable and portable.

Dimensions: L 7.5" x W 4" x D 3"

Weight: 1 lb.

Technology

The Mobilize Compact bleeding control and airway management system includes a multi-user license for the Mobilize Rescue app, which users access on their smartphone (with or without WiFi). The app provides real-time treatment instructions through interactive prompts, written instructions, and visual cues. It is available for both iOS and Android devices.



The Mobilize Rescue app provides real-time treatment guidance and is available for iOS and Android devices.

Equipment and App Integration

All medical supplies are clearly labeled, logically organized, and color-coded to match the instructions in the Mobilize Rescue app. The system's intuitive organization and integrated instructions makes locating supplies easy, saving valuable time during an emergency.



Training Kits

The Mobilize Training Kit includes the same supplies and app as the Compact unit. Use it to familiarize people in your organization with the rescue systems and app functionality and benefits.

STORAGE RECOMMENDATIONS & DIAGNOSTIC CHECK | Individual supplies within the unit may be negatively impacted by exposure to hot and cold temperatures. Individual supplies' expiration dates vary from 1-5 years from date of purchase; it is crucial that you perform diagnostic checks to ensure that the unit is charged to 100%, the proof seal is intact (not tampered with), and that the medical supplies are not damaged or contaminated.

¹Mell, H., Mumma, S., Hiestand, B., Carr, B., Holland, T., Stopyra, J. (2017, July 19). Emergency medical response times in rural, suburban and urban areas [Editorial]. JAMA Surgery DOI: 10.1001/jamasurg.2017.2230

Copyright ©2021 ZOLL Medical Corporation. All rights reserved. Mobilize and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

MCN PP 2107 0598

ZOLL Medical Corporation
Worldwide Headquarters
269 Mill Road
Chelmsford, MA 01824, USA
978-421-9655
800-804-4356

For subsidiary addresses and fax numbers as well as other global locations, please go to zoll.com/contacts.

ZOLL