

ZOLL Intravascular Temperature Management Catheters

ZOLL's portfolio of catheters provides different lengths, insertion sites, and heat-exchange power to meet patient needs, clinical scenarios, and protocols. ZOLL offers four intravascular temperature management catheter choices that connect to the Thermogard XP® console.

All Catheters Feature:

- **Triple-lumen central venous catheter (CVC) functions**, including medication delivery, blood draw, and central venous pressure monitoring
- Accessories needed for placement
- Hydrophilic coating with heparin
- Radiopaque body, tip, and marker band to ensure proper placement in the vessel

All Catheters:

- 9.3F
- MR conditional
- Do not contain natural rubber latex

SOLEX 7™ CATHETER

- For fever management and perioperative temperature management, with double the power of Cool Line® (when used with Thermogard XP®)
- 7-day dwell time for fever management and 4 days for perioperative temperature management
- Internal jugular or subclavin insertion sites
- Serpentine balloon design
- 20 cm length



Description

Catalog Number

Solex 7 Catheter Kit

8700-0793-01

Solex 7 Disposable Pack: Includes four (4) Solex 7 Catheter Kits and four (4) Start-Up Kits (6 ft)

8700-000868-01



QUATTRO® CATHETER

- For reaching/maintaining normothermia during surgery and recovery/intensive care
- 4-day dwell time
- Femoral insertion site
- 4 heat-exchange balloons
- 45 cm length

Description

Catalog Number

Quattro Catheter Kit

8700-0783-01

Quattro Disposable Pack: Includes four (4) Quattro Catheter Kits and four (4) Start-up Kits (6 ft)

8700-000867-01

Quattro Catheter Kit, with Premium Access Kit

8700-0791-01

Quattro Premium Disposable Pack: Includes four (4) Quattro Catheter Kits and four (4) Start-up Kits. Each catheter kit contains insertion accessories including full body drape, ChloroPrep® applicator, and retractable scalpel.

8700-000766-01





COOL LINE® CATHETER

- For fever management
- 7-day dwell time
- Subclavian, internal jugular, or femoral insertion sites
- 2 heat-exchange balloons
- 22 cm length

Description

Cool Line Catheter Kit

Catalog Number

8700-0781-01

Cool Line Disposable Pack: Includes four (4) Cool Line Catheter Kits and four (4) Start-Up Kits (6 ft).

8700-000865-01



ICY® CATHETER

- For reaching/maintaining normothermia during surgery and recovery/intensive care
- 4-day dwell time
- Femoral insertion site
- 3 heat-exchange balloons
- 38 cm length



Description

Icy Catheter Kit

Catalog Number

8700-0782-01

Icy Disposable Pack: Includes four (4) Icy Catheter Kits and four (4) Start-up Kits (6 ft)

8700-000866-01



START-UP KIT

- Start-up Kit, standard tubing: 6 feet
- Start-up Kit, extended tubing: 9 feet



Description

Start-up Kit, standard tubing (6 feet)
Start-up Kit, extended tubing (9 feet)

Catalog Number

8700-0784-01
8700-0785-01

To learn more about these products, contact your local ZOLL Temperature Management representative or the ZOLL Customer Service Department at 1-800-348-9011.