

ELECTRODE REFERENCE GUIDE

Product Name & Number	Differentiation Features	Primarily Used In	Defib Compatibility		
 stat•padz (8900-4003)	<ul style="list-style-type: none"> Green pressure sponge - Promotes better skin coupling. Pre-Connect - Saves critical application time. Speed pack - Organizes wires to eliminate tangles and fumbling. Different pad shapes - Visual cues for proper pad placement. 	<ul style="list-style-type: none"> Med/surge floors ER & ICU EMS & Fire 			
 Radiolucent Solid Gel (8900-4005)	<ul style="list-style-type: none"> For use in Cath and EP labs for procedures - Virtually invisible to X-ray. Pre-Connect - Saves critical application time. Speed pack - Organizes wires to eliminate tangles and fumbling. Different pad shapes - Visual cues for proper pad placement. 	<ul style="list-style-type: none"> Cardiac catheterization lab Electrophysiology labs ER & ICU 			
 Radiolucent Liquid Gel (8900-2105-01)	<ul style="list-style-type: none"> For use in Cath and EP labs for procedure. Liquid gel provides optimal skin coupling for no dermal damage. Virtually invisible to X-ray. 	<ul style="list-style-type: none"> Cardiac catheterization lab Electrophysiology labs 			
 Sterile pro•padz (8900-4012)	<ul style="list-style-type: none"> For use in operating rooms. Helps maintain sterile field. Packaged in two pouches. 	<ul style="list-style-type: none"> Operating rooms 			
 Cardiology Specialty (8900-2100-01)	<ul style="list-style-type: none"> No-slump tray. Liquid gel provides optimal skin coupling for no dermal damage. 	<ul style="list-style-type: none"> Cardiac catheterization lab Electrophysiology labs 			
 Biphasic pro•padz (8900-2303-01)	<ul style="list-style-type: none"> Used in Cath and EP labs for procedures. Highest conversion rate for cardioversions with low energy biphasic waveform. Liquid gel provides optimal skin coupling for no dermal damage. 	<ul style="list-style-type: none"> Cardiac catheterization labs Electrophysiology labs 			
 pedi•padz (8900-2065)	<ul style="list-style-type: none"> Meets all pediatric ACLS criteria. Liquid gel provides optimal skin coupling for no dermal damage. 	<ul style="list-style-type: none"> Med/surge floors ER & ICU EMS & Fire 			
 pedi•padz Solid Gel (8900-3000-01)	<ul style="list-style-type: none"> Meets all pediatric ACLS criteria. Solid gel provides ease for storing. Pre-connect saves critical application time. 	<ul style="list-style-type: none"> EMS & Fire 			
 pedi•padz Radiolucent (8900-1005-01)	<ul style="list-style-type: none"> Meets all pediatric ACLS criteria. For use during procedures. Virtually invisible to X-ray. 	<ul style="list-style-type: none"> Cardiac catheterization lab Electrophysiology labs 			
 pedi•padz II (8900-0810-01)	<ul style="list-style-type: none"> Gives AEDs Pediatric intelligence function. Pediatric intelligence function only available on ZOLL® AEDs with pedi•padz® II. 	<ul style="list-style-type: none"> EMS & Fire Public areas (i.e. lobby) 			 
 stat•padz II (8900-0802-01)	<ul style="list-style-type: none"> Connects directly to AEDs and ALS devices. See stat•padz for other benefits. 	<ul style="list-style-type: none"> EMS & Fire 			 
 CPR-D•padz (8900-0800-01)	<ul style="list-style-type: none"> Provides defibrillator with Real CPR Help®. One-piece design simplifies training and retention of training, and assures proper positioning on the chest. One size fits all, with ability to expand for larger chest. 	<ul style="list-style-type: none"> EMS & Fire Public Areas (i.e. lobby) 	See Reverse for more info	See Reverse for more info	See Reverse for more info 
 CPR stat•padz (8900-0400)	<ul style="list-style-type: none"> Designed for ALS market, two-piece design ensures all therapies can be delivered. Provides defibrillator with Real CPR Help. Accommodates AP and AA pad placement. 	<ul style="list-style-type: none"> Med/surge floors ER & ICU EMS & Fire 	See Reverse for more info	See Reverse for more info	See Reverse for more info 
 OneStep Electrodes	<ul style="list-style-type: none"> OneStep™ Basic 8900-0211-01 OneStep Pacing 8900-0212-01 OneStep CPR 8900-0213-01 One Step CPR (AA) 8900-0217-01 OneStep Complete 8900-0214-01 <p>* See OneStep electrode spec sheet for more information</p>	<ul style="list-style-type: none"> Med/surge floors ER & ICU 	See Reverse for more info	See Reverse for more info	

ZOLL CPR ELECTRODE SOLUTIONS


CPR-D•padz
(8900-0800)

CPR stat•padz Case
(8900-0400)

OneStep CPR
(8900-0213-01)

OneStep CPR AA
(8900-0217-01)

OneStep Complete
(8900-0214-01)



AEC PLUS



- Connects directly for full CPR capabilities

- Connects directly for full CPR capabilities

- Does not connect

- Does not connect

- Does not connect

AEDPRO



- Connects directly for full CPR capabilities

- Connects directly for full CPR capabilities

- Does not connect

- Does not connect

- Does not connect

M Series



- Needs CPR Connector (pt# 8000-0370) for CPR
- Connects directly for defib

- Needs CPR Connector (pt# 8000-0370) for CPR
- Connects directly for defib/pace

- Needs new M Series® CPR Cable for CPR (pt# Not Available)
- Connects directly for defib/pace

- Needs new M Series CPR Cable for CPR (pt# Not Available)
- Connects directly for defib/pace

- Needs new M Series CPR Cable for CPR (pt# Not Available)
- Connects directly for defib/pace

E SERIES



- Needs CPR Connector (pt# 8000-0370) for CPR
- Connects directly for defib

- Needs CPR Connector (pt# 8000-0370) for CPR
- Connects directly for defib/pace

- Does not work for CPR
- Connects directly for defib/pace

- Does not work for CPR
- Connects directly for defib/pace

- Does not work for CPR
- Connects directly for defib/pace

R series™



- Needs R Series® CPR Connector for CPR (pt# Not Available)
- Connects directly for defib

- Needs R Series CPR Connector for CPR (pt# Not Available)
- Connects directly for defib/pace

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

Function	Defibrillate Only	Multifunction	Multifunction	Multifunction	Multifunction
Construction	1 piece	2 piece	2 piece	2 piece	2 piece
Pad Placement	AA	AA/AP	AP	AA	AP
Shelf Life	5 yrs	2 yrs	2 yrs	2 yrs	2 yrs
Quantity Sold	Individual	8 per Case	8 per Case	8 per Case	8 per Case