AED Plus[®]Accessories and Training Products





ACCESSORIES FOR EVERY RESCUE

ZOLL ResQSystem[™] Cabinet

One location. Many ways to help.

Part # 8911-005050

Spacious, centrally located cabinet for AEDs and trauma kits

- Quick and easy access
- AED and medical-supply graphic indicators on front and sides
- Charging port for Mobilize Comprehensive trauma kit ensures device is charged
- Suitable for all types of organizations





Semi-recessed Wall Cabinet

Part # 8000-001257

Meets ADA (Americans with Disabilities Act) requirements

Dimensions: $14'' \times 14'' \times 3.5''$



Fully Recessed Wall Cabinet

Part # 8000-001258

Meets ADA requirements

Dimensions: $14'' \times 14'' \times 1.5''$



9-inch Standard Wall Cabinet

Part # 8000-0855

Alarmed wall cabinet holds the AED Plus defibrillator with the carry case and one spare set of electrodes.

Dimensions: $17.5" \times 17.5" \times 9"$



7-inch Standard Wall Cabinet

Part # 8000-001256

Alarmed wall cabinet holds the AED Plus defibrillator *without* the carry case and one spare set of electrodes.

Dimensions: $17.5" \times 17.5" \times 7"$



Strobe for Standard Surfacemounted Wall Cabinet

Part # 8000-001259



Strobe for Recessed Wall Cabinet

Part # 8000-001267

Electrodes, Batteries & Accessory Kit

Everything you need to keep your AED Plus ready to use.



CPR-D-padz® Electrodes

Part # 8900-0800-01

- One-piece design ensures fast and accurate placement
- For use on adult patients and children over the age of 8 or > 55 pounds (25 kilograms)
- Built into the electrodes, Real CPR Help® technology senses the rate and depth of chest compressions
- Five-year shelf life



Accessory Kit

Part # 8900-0807-01

• The AED Plus accessory kit includes scissors, disposable razor, face barrier/shield, and dry towel.

Note: This kit is included and attached to every set of ZOLL CPR-D-padz® electrodes.



Pedi-padz® II Pediatric Electrodes

Part # 8900-0810-01

- ZOLL's intelligent pediatric capability is built into the AED Plus defibrillator. The AED can sense when Pedi-padz II are being used and will switch from adult to pediatric mode, modifying the device to provide an appropriate pediatric ECG analysis algorithm and therapy energy level.
- For use with patients under the age of 8 years or < 55 pounds (25 kilograms)
- Two-year shelf life

Caution: Pediatric electrodes should never be used with an adult patient.



Type 123 Lithium Batteries (sleeve of 10)

Part # 8000-0807-01

- Non-proprietary off-the-shelf Duracell
 123 lithium batteries
- 5-year standby life once installed in AED Plus

Caution: Only use batteries from a ZOLL-approved manufacturer (refer to User Manual).

Mounting Bracket & Signage

Ensure AED visibility



Mounting Bracket
Part # 8000-0809-01

Bracket for wall-mounting AED Plus without carry case



AED Wall Decal

Part # 8000-0849-01

4" x 4" AED decal suitable for indoor or outdoor use



AED Sign Package

Part # 8000-0825

3D wall sign and 8.5" x 11" wall-mounted sign to improve AED visibility

Carry Cases, Backpacks, & Bags

Options for storage and transport



Rescue Backpack

Part # 8000-001468-01

- Outer pouch that holds the AED Plus
- Three removeable, stick-in-place pouches in main compartment
- Side pockets for easy access to smaller equipment

Dimensions: 19" x 18" x 8"



Large Hard-shell Carry Case

Part # 8000-0837-01

- Protection for AED Plus in more demanding environments
- Large enough to store AED Plus, CPR-D-padz, Pedi-padz
 II, and sleeve of 10 batteries

Dimensions: $19'' \times 15.4'' \times 7.6''$



Small Hard-shell Carry Case

Part # 8000-0836-01

- Protection for AED Plus in more demanding environments
- Stores AED Plus defibrillator only

Dimensions: 16" x 13" x 6.9"



Soft Carry Bag Part # 8000-0802-01

- Provides easy access and quick retrieval of AED Plus and spare accessories
- Pockets for storing extra electrodes and batteries



AED Device Management & Upgrade Tools

Administrative Software (ZAS)

Part # 8000-0311-01

- Software to reconfigure settings on ZOLL AED Plus
- Requires use of IrDA adapter (sold separately)
- Software can be downloaded at no cost directly from ZOLL website

AHA Guidelines Upgrade Kit

Part # 7771-000011-01

- Upgrades ZOLL devices manufactured prior to 2010 to current AHA Guidelines
- Includes Administrative Software (ZAS) CD, instructions, and labels to upgrade AED Plus faceplate
- Requires use of IrDA adapter (sold separately)

USB-to-IrDA Adapter Kit

Part # 8000-0815

- USB-to-IrDA adapter cable and driver software
- Required for AED Plus configuration and post-rescue event data review

Mini USB-to-IrDA Adapter Kit

Part # 8000-1200

- Mini USB-to-IrDA adapter cable and driver software
- Required for AED Plus configuration and post-rescue event data review

RS232-to-IrDA Adapter Kit

Part # 8000-0816

- RS232-to IrDA adapter cable and driver software
- Required for AED Plus configuration and post-rescue event data review

RescueNet® EventSummary Software

- Provides a report of ECG and vital sign data, CPR rate and compression depth, and all device actions and voice prompts
- Summary reports exportable as PDFs

- Detailed information allows for complete review and analysis of the rescue
- Easy to email to patient's care provider or print for post-event debriefing and education

Download for free at: www.zoll.com/rescueneteventsummary



AED Plus Clinical Unit Training

Devices and Accessories

AED Plus Demo Kit

Part # 8000-0834-01

- Kit to demonstrate Real CPR Help® with clinical AED Plus
- Includes CPR-D Demo electrodes and manikin in large carry bag
- Clinical AED Plus and simulator not included





ZOLL AED Simulator Part # 8000-000925

 Simulates heart rhythm for AED Plus demonstration using clinical unit



Travel Trainer

Part # 8008-0006-01

AECPLUS

- One-piece, compact, easy-to-carry demonstration tool
- Simulates VF and normal heart rhythms



Demo Manikin*

Part # 8000-0835-01

- Adult manikin for AED/CPR training
- Use with AED Plus clinical unit and AED Plus Trainer2
- Training electrodes and simulator sold separately

AED Plus Trainer 2 & Accessories

Everything you need to simulate a rescue with your AED Plus Trainer2

AED Plus Trainer2

Semiautomatic: Part # 8008-0050-01 Fully automatic: Part # 8008-000052-01

- Allows the student to experience what it's like to help save a life with the AED Plus defibrillator. Choose from four different training scenarios or have complete manual control of the entire rescue. A wireless remote lets you control multiple Trainer2 devices.
- Ships with Operator's Manual, batteries, wireless remote, and adult training electrodes
- Replacement adult training electrodes, pediatric training electrodes, and adhesive gel can be purchased separately as needed





Pedi-padz II Training Electrodes

Part # 8900-000861-01

- Pediatric electrodes for use with the Trainer2
- Pack of 6



Demo Manikin*

Part # 8000-0835-01

- Adult manikin for AED/CPR training
- Use with AED Plus clinical unit and AED Plus Trainer2
- Training electrodes and simulator sold separately



CPR-D-padz Training Electrodes

Part # 8900-0804-01

 Adult electrodes for use with the Trainer2. Replacement gel can be purchased separately.



Replacement Adhesive Gel for CPR-D-padz Training Electrodes

Part # 8900-0803-01

- Gel must be applied to the reusable pad before use.
- Pack of 5



AED Plus Trainer2 Soft Carry Case

Part # 8000-0375-01

 Clearly identifies the device as a training unit and provides a convenient way to transport the Trainer 2. Includes a pocket for accessories



Replacement Trainer2 Wireless Remote

Part # 8008-0007

 Wireless remote for use with the Trainer2



AED Program Management

Stay prepared to save a life

PlusTrac® Professional AED Program Management for a single AED

1 Year: Part # 8000-1110-01 5 Years: Part # 8000-1111-01

- Medical authorization/Rx
- Online management software access
- Protocol and EMS templates
- AED law library and support knowledge base
- Email reminders and alert notifications
- Complete customer support/implementation
- Post-event services and loaner program

PlusTrac Essentials AED Program Management for a single AED

1 Year: Part # 8000-1116-01 5 Years: Part # 8000-1117-01

- Medical authorization/Rx
- Online management software access
- Protocol and EMS templates
- AED law library and support knowledge base
- Email reminders and alert notifications



For general product-related questions or assistance purchasing ZOLL accessories, contact Customer Support at 800-348-9011. For technical questions about your ZOLL AED Plus defibrillator or accessories, contact Technical Support at 978-421-9460.

ZOLL MEDICAL CORPORATION

269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | 800-804-4356 | zoll.com

For subsidiary addresses and fax numbers, as well as other global locations, please go to zoll.com/contacts.



Copyright © 2024 ZOLL Medical Corporation. All rights reserved. AED Plus, CPR-D-padz, Pedi-padz, PlusTrac, Real CPR Help, RescueNet, ResQSystem, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners. Printed in U.S.A.